

Rapicide® OPA/28

HIGH-LEVEL DISINFECTANT



The Fastest and Longest Lasting OPA on the Market Today!

The RAPICIDE OPA/28 Advantage **BY THE NUMBERS**

- 5** MINUTE AER disinfection time at 25°C
- 10** MINUTE MANUAL disinfection time at 20°C minimum
- 84** Up to 10-24 additional AER reuse cycles†
- 28** DAY reuse period, twice as long as other OPA brands
- 2** Final Rinses needed in MEDIVATORS® AER
- 3** Final Rinses needed in manual applications



WORKS FASTER

Take quicker action against microorganisms



High-level disinfects devices manually in just 10 minutes at room temperature (20°C minimum), to achieve bactericidal, fungicidal and tuberculocidal kills. Studies indicate Rapicide OPA/28 effectively inactivates:

MRSA • VRE • TB
Hepatitis Viruses • HIV

OPA's faster manual disinfection allows you to reprocess one extra scope per hour as compared to other OPA products.

LASTS LONGER

Use for up to 28 days

It's important for healthcare providers to manage budgets carefully while still delivering high-quality patient care. Rapicide OPA/28 has been tested to provide 84 reuse cycles† along with a reuse period of up to 28 days – twice the reuse period of other OPA brands. It also features a 2-year unopened shelf life and 75 days once opened.

† MEDIVATORS® International Testing Per EPA Reuse Test Protocol Specification.



Distributed by:



www.supply4gi.com

Description	Quantity	Cat. No.
Rapicide® OPA 28 Day HLD	1 Gallon Bottle	S4GIRapicideOPA
Also Available: Rapicide® OPA Test Strips	50 Strip Bottle	S4GIROPATest